

# Amvisc®

1.2% Sodium Hyaluronate



VISCOELASTIC

A general purpose viscoelastic that enables excellent visualization and chamber maintenance throughout the procedure.

# OPTICAL

CLARITY

# PRESERVE **TISSUE** INTEGRITY



- Creates and maintains a deep chamber for a controlled capsulorhexis and clear, unobstructed view of the rhexis perimeter.<sup>1</sup>

- Preserves tissue integrity by shielding corneal endothelial cells from surgical stresses.<sup>1</sup>

**Amvisc®**  
1.2% Sodium Hyaluronate

**BAUSCH + LOMB**

# Amvisc®

## 1.2% Sodium Hyaluronate


LEARN MORE AT [bauschsurgical.com](http://bauschsurgical.com)



59081L order number 59081LXXXX



MODEL NUMBER	SIZE	VISCOSITY, CP (AT 1S-1 SHEAR RATE, 25°C)	MOLECULAR WEIGHT (DALTONS)	COMPOSITION	COHESION	OSMOLALITY (MOSM/KG)	CANNULA SIZE
59081L	1.0mL	40,000	>1,000,000	1.2%HA	Cohesive	320	27G



### PRODUCT NUMBER 59081L

#### DESCRIPTION

- High-viscosity, 1.2% sodium hyaluronate (HA) concentration
- 0.8-mL preloaded sterile glass syringe with finger grip for comfort
- 27-gauge disposable cannula, 32° angle
- Syringe requires assembly of the cannula only
- Not made with natural rubber latex
- >1.0 million molecular weight (Daltons)

#### DISPENSED

Singles

1. Amvisc® Directions for Use



Find B+L IOL surgical equipment online at [www.StorzEye.com](http://www.StorzEye.com)

## Indications and Important Safety Information for Amvisc® and Amvisc® Plus Viscoelastics

**INDICATIONS:** Amvisc® and Amvisc® Plus are indicated for use as a surgical aid in ophthalmic anterior and posterior segment procedures including: Extraction of a cataract, implantation of an intraocular lens (IOL), corneal transplantation surgery, glaucoma filtering surgery, and surgical procedures to reattach the retina.

Because of its lubricating and viscoelastic properties, transparency, and ability to protect corneal endothelial cells, Amvisc and Amvisc Plus helps maintain anterior chamber depth and visibility, minimizes interaction between tissues, and acts as a tamponade and vitreous substitute during retinal reattachment surgery. Amvisc and Amvisc Plus also preserves tissue integrity and good visibility when used to fill the anterior and posterior segments of the eye following open sky procedures.

**PRECAUTIONS:** There may be increased intraocular pressure following surgery caused by pre-existing glaucoma or by the surgery itself. An excess quantity of Amvisc and Amvisc Plus should not be used. Amvisc and Amvisc Plus should be removed from the anterior chamber at the end of surgery to prevent or minimize post-operative intraocular pressure increases (spikes). If the postoperative pressure increases above expected values, correcting therapy should be administered. Reuse of cannula should be avoided. Store at 2-8° C. Protect from freezing.

**ADVERSE EVENTS:** Transient postoperative inflammatory reactions were reported in clinical trials and oral and topical steroid preparations were administered. The best index of the degree of phlogistic response is the postoperative clarity of the vitreous cavity. Transient postoperative increase in intraocular pressure has been observed following the use of sodium hyaluronate in anterior segment surgery. On rare occasions postoperative reactions, including inflammation, corneal edema, and corneal decompensation have been reported.

**ATTENTION:** Refer to the Directions for Use labeling for a complete listing of indications, warnings and precautions, clinical trial information, etc.

**CAUTION:** Federal (USA) law restricts this device to the sale by or on the order of a physician.