

Single use instrument, including handle  
 No cleaning or sterilization required  
**COMPLETELY DISPOSABLE**

CapsuleGuard<sup>®</sup> IA

Irrigation/Aspiration Handpiece






EXCELLENT  
 FOR USE IN  
**ALL**  
 PHASES OF I/A

AVAILABLE IN 3 TIP STYLES & ADAPTABLE WITHIN A LARGE RANGE OF INCISION SIZES


Stellaris<sup>®</sup> System



		MICS 1.8mm	Standard Incision 2.2-2.8mm
Straight		85914ST	85912ST
Curved			85913ST
45°		85915ST	85910ST

Other Systems



		MICS 1.8mm	Standard Incision 2.2-2.8mm
Straight		85914S	85912S
Curved			85913S
45°		85915S	85910S

PRODUCT FEATURES

- Flexible Tip
- Two Irrigation Ports
- One Piece Silicone Construction
- Flexible Irrigation Tube
- Soft Aspiration Port
- Flexible Sleeve Conforms to Wound
- Enhanced Visualization from Semi-transparent Silicone
- No Sharp Edges for Ease of Insertion Into the Eye

**STORZ**<sup>®</sup>  
 Ophthalmic Instruments

800-338-2020  
[www.StorzEye.com](http://www.StorzEye.com)