

DISPOSABLE


# Cohen-Benner Contact Lens

Flat contact lens for  
posterior visualization



Designed to grip the edge of the cornea, offering a free-floating, 2-piece design for contact with the tear film.

**FREE-  
FLOATING**



LENS DESIGN THAT ADJUSTS TO  
DIFFERENT CORNEAL TOPOGRAPHIES

**SMALL  
FOOTPRINT**  
WITHOUT COMPROMISING 36 DEGREE FOV

PROPRIETARY  
**STABILIZING** SYSTEM

4 STUDS FOR CAPILLARY TRACTION



Exclusively available from **BAUSCH+LOMB**

# Cohen-Benner Contact Lens

MODEL#54013



A set-it-and-forget-it contact lens



## DDMS™ Tano

- 20.04 20G Tano
- 20.04.23 23G Tano
- 20.04.25S 25G Tano
- 20.04.27 27G Tano

☐=5 STERILE ⓘ



## DDMS™ Retractable Tano

- 20.16.23 23G Retractable Tano
- 20.16.25 25G Retractable Tano

☐=5 STERILE ⓘ



## Fine Tip Eckardt Forceps

- 20G 09.10
- 23G 0.9.10.23PIN; D39.10.23; 09.10.23
- 25G 0.9.10.25PIN; D39.10.25; 09.10.25

☐=5 STERILE ⓘ



## Microserrated™ Tewari Omnicap Forceps

- 23G 39.15.23PIN, D39.15.23
- 25G 39.15.25PIN, D39.15.25
- 27G 39.15.23PIN

☐=5 STERILE ⓘ



## Soft Tip Cannula

- 30.08.23 23G
- 30.08.25 25G
- 30.08.27 27G

☐=6 STERILE ⓘ



**CONTACT YOUR B+L REPRESENTATIVE  
TO PLACE YOUR ORDER**

Click to learn more.