

freeflow™ infusion

The unique cap-over configuration of the FreeFlow™ infusion line is designed to complement the efficiency of the Stellaris Elite™ platform and Bi-Blade® vitrectomy cutter. It's distinctive silicone sleeve snaps over the head of the cannula, allowing a higher volume of fluid to pass through.* This facilitates a stable environment, even at higher flow rates and cutting speeds up to 15,000 CPM.


Additionally, its unique 30° angled design drapes cleanly along the contour of the patient's face compared to traditional infusion lines, keeping more of your workspace obstruction free.

*Compared to standard infusion lines utilizing a cannula inside a cannula design.

Ordering Information

Product #	Description (6/Box)
30.35	FreeFlow Infusion

Flow Rates*
Standard vs. FreeFlow Infusion Lines



*Average water flow (mL/min) at an infusion pressure of 30mmHg

Engineered for

Stellaris Elite
Vision Enhancement System

bi-Blade

Please see www.bauschretina.com for Important Safety Information.

BAUSCH+LOMB



Synergetics