



I-PACK™

Injection Kit



Choose Convenience for
Your Intravitreal Injections

ALL-IN-ONE

COMPLETE & STERILIZED KIT

STREAMLINE THE PROCESS

OF PERFORMING INTRAVITREAL INJECTIONS

At Bausch + Lomb, we strive to provide you with innovative products that help meet your ever-evolving needs. I-PACK™ delivers a convenient solution that streamlines the process of performing intravitreal injections. Packaged in a complete and sterilized kit, this all-in-one solution has everything you need, when you need it. It can help you reduce ordering time, organize inventory and simplify preparation time.

I-PACK™ components are delivered within the same plastic tray that serves as the reservoir for your povidone-iodine solution. Simply pour the solution into the compartment where it can be used for dipping cotton swabs. To make cleanup one simple step, the two additional compartments can be used for discarding vials, swabs and needles.

BAUSCH + LOMB

I-PACK™ Injection Kit



STREAMLINE ORDERING
ORGANIZE INVENTORY
SIMPLIFY PREPARATION

I-PACK™ Injection Kit Custom comes in a variety of pack configurations so you can customize for your unique needs.

18060 Components

- (1) 1 mL Tuberculin Syringe
- (1) 30ga 1/2" Needle
- (1) 27ga 1/2" Needle
- (2) 18ga 1" Needles
- (6) Cotton Swabs 6"
- (1) Barraquer Eyelid Speculum, Adult
- (1) 3.5-4mm Scleral Marker

18061 Components

- (2) 1 mL Tuberculin Syringes
- (1) 30ga 1/2" Needle
- (1) 27ga 1/2" Needle
- (2) 18ga 1" Needles
- (6) Cotton Swabs 6"
- (1) Barraquer Eyelid Speculum, Adult
- (1) 3.5-4mm Scleral Marker
- (1) Eye Drape (16" x 16")

18064 Components

- (2) 1 mL Tuberculin Syringes
- (1) 30ga 1/2" Needle
- (1) 27ga 1/2" Needle
- (2) 18ga 1" Needles
- (6) Cotton Swabs 6"
- (1) Barraquer Eyelid Speculum, Adult
- (1) 3.5-4mm Scleral Marker
- (1) Eye Drape (16" x 16")
- (1) Disposable Forceps, 4"
- (4) 4" x 4" Gauze Pad

18067 Components

- (1) 30ga 1/2" Needle
- (1) Barraquer Eyelid Speculum, Adult
- (1) 3.5-4mm Scleral Marker

18069 Components

- (1) 30ga 1/2" Needle
- (6) Cotton Swabs 6"
- (1) Barraquer Eyelid Speculum, Adult
- (1) 3.5-4mm Scleral Marker

18069W Components

- (1) 30ga 1/2" Needle
- (6) Cotton Swabs 6"
- (1) Wire Eyelid Speculum, Adult
- (1) 3.5-4mm Scleral Marker



CONTACT YOUR B+L REPRESENTATIVE TO PLACE YOUR ORDER

Click to learn more.