

# FOREIGN BODY



Specialized design for unique foreign body cases with a glare reducing finish for better visualization.

## BB Removal Forceps

18G 11.40PIN =5 STERILE

- Designed for the removal of a 4.5mm ball bearing (BB)



## Foreign Body Forceps

19G 11.20PIN =5 STERILE

- Carbide coated grasping surface with 2.5mm tip opening



## IO Foreign Body Magnet

18G 12.21PIN =5 STERILE

- Capable of lifting a 11/16" (17.5mm) steel ball bearing
- Actuate instrument to shield magnet and release foreign body



Engineered for



Please see [bauschurgical.com](http://bauschurgical.com) for Important Safety Information.

**BAUSCH+LOMB**



**Synergetics®**