

REDDY END GRASPING FORCEPS



External surfaces of the distal tip are micro serrated to assist in edge identification and initiate membrane peels

Reddy End Grasping Forceps

23G 39.32.23PIN =5  STERILE 

25G 39.32.25PIN =5  STERILE 

- Tapered tip provides a precise grasp
- Long, micro serrated grasping platform to maintain purchase of ILM/ERM
- Glare reducing finish for better visualization



Engineered for



Please see bauschsurgical.com for Important Safety Information.

BAUSCH+LOMB



Synergetics®